

Human Serum AB male HIV tested Heat Inactivated

Collected from the source:

When searchers choose their serum, an important factor that should be taken into consideration is the source, which also emphasises the traceability of the serum.

Our system of vertical integration allows us to be certain of the origins and traceability of our human male AB serum. The donors are volunteers.

Each manufactured batch is rigorously controlled, from the collection of serum in authorised organisms, and throughout all stages of its treatment and production through to final packaging on the authorised organisms and on our premises.

The serum is off-the-clot serum, processed from human male AB group blood that has had natural coagulation.

The serum is collected or imported and treated in agreement with the European regulations.

Filtration:

Final Filter Size: 0.2µm x 2

Sterility:

All sera are tested for the absence of aerobic and anaerobic bacteria, fungi and yeast.

The sterility test is based on the European Pharmacopeia requirements.

The sera are tested for the absence of *Mycoplasma pneumoniae* by chimiluminescence.

Virus and antibodies tested:

All of our human serum is tested for:

- Hepatitis B antigen (HBs Ag)
- Hepatitis C virus and antibodies (HCV)
- HIV Type 1 virus and antibodies HIV ½
- Syphilis

Endotoxin:

All sera are tested to determine the levels of endotoxins. Serox carries out a chromokinetic quantitative test, according to the method D of the European Pharmacopoeia.

The endotoxin reagent is standardized against the US reference endotoxin.

Haemoglobin:

The haemoglobin level is measured by spectrophotometer.

Osmolality:

Determined by a lowered freezing temperature. The osmometer is calibrated against standard solutions.

Cell Culture:

Biological performance is assessed using cell culture medium supplemented with the serum being tested.

During the test period, cultures are examined microscopically for any morphological abnormalities that may indicate toxic components in the serum.

Cell Lines Tested:

MRC-5 - Human Fetal Lung

HELA - Cancer Cell/Human

Total Protein:

Determined by Biuret Colorimetry.

Country of Origin:

It is the country in which the serum was taken from the donor.
Serox sera are sourced from France, Germany, Poland or USA.

Storage conditions:

Store at -20°C.

Shelf life:

5 years.

Recommended use:

- Respect storage conditions of the serum
- Do not use the serum after its expiry date
- Store serum in an area protected from light
- Manipulate serum in aseptic conditions (e.g.: under laminar air flow)
- Wear clothes adapted to the manipulation of serum to avoid contamination (e.g.: gloves, mask, hygiene cap, overall...)
- In order to preserve all serum qualities, it is recommended to thaw out the flask, to aliquote, then to re-freeze the produced flasks rather than to thaw out and re-freeze the flask at each use.
- It is recommended to use the serum immediately after its thaw out. However, if it is not useful, it is possible to store thaw out serum, at +2°C / +8°C, until 26 weeks without significant decrease of its performances in cell culture.

The product is intended to be used in vitro, in laboratory only. Do not use it in therapy, human or veterinary applications.

Heat inactivation :

Sera are heat inactivated to inactivate complement. The complement can lead to complement-mediated cell lysis. Immunological studies justify the need to heat inactivate the serum. The treatment is a heating of the serum at 56°C during 30 minutes.

This treatment can modify the colour of the serum. This is normal and it does not compromise the quality of the product for the cell culture.

Note : To heat the serum during a long period can reduce or destroy the growth factors. It can also increase the build of precipitates that are frequently confounded with contamination.

Effects of the heat inactivation :

- Destruction of proteins
- Precipitation of fibrin
- Destruction of fibrinogen
- Deterioration total or partial of vitamins
- Decrease of the growth factors concentration
- Destruction of the LDH
- Decrease of the amylase concentration
- Destruction of the alkaline phosphatase
- Deterioration of the IgG, IgE, IgM
- Inactivation of viruses
- Inactivation of T4 Phage
- Inactivation of mycoplasma
- Increase of oxidation and catalysis reactions.

The chemical composition of the serum is not altered.