

Ref: FT.L0191an Version date: 29/02/16

Iscove's Modified Dulbecco's Medium

w/ Stable Glutamine w/ 25mM Hepes

CAT N°: SLIM-537-500

Theoretical pH : 7.2 ± 0.3

Osmolality: 290 mOsm/kg \pm 10 %

Colour: Red, clear solution

Storage conditions : $+2^{\circ}$ C to $+8^{\circ}$ C. Protect from light.

Shelf life: 24 months

Sterility tests:

- Bacteria in aerobic and anaerobic conditions

- Fungi and yeasts

Endotoxin: < 1 EU/ml

Cell growth test:

Medium tested for the ability to support SP2/0-Ag14 cell growth.

Composition: Displayed on website and in catalogue; also available on request.

Recommended use:

- Respect storage conditions of the product
- Do not use the product after its expiry date
- Store product in an area protected from light (not necessary for saline solutions).
- Manipulate the product in aseptic conditions (e.g. : under laminar air flow)
- Wear clothes adapted to the manipulation of the product to avoid contamination (e.g. : gloves, mask, hygiene cap, overall...)

The product is intended to be used in vitro, in laboratory only. Do not use it in therapy, human or veterinary applications.

Application:

In 1976, Guilbert and Iscove demonstrated that precursor cells of erythrocytes and macrophages could be cultured in a reduced-serum medium supplemented with albumin, transferrin, lecithin, and selenium.

Iscove's medium is a modification of Dulbecco's Modified Eagle's Medium (DMEM) containing selenium, additional amino acids and vitamins, sodium pyruvate, HEPES buffer, and potassium nitrate instead of ferric nitrate .

Further studies demonstrated that Iscove's Medium would support murine B lymphocytes, hemopoietic tissue from bone marrow, B cells stimulated with lipopolysaccharide, T lymphocytes, and a variety of hybrid cells.

The stable Glutamine prevents the intramolecular cyclization reaction associated with solutions of L-Glutamine. This derivative is therefore stable in solution and allows the formulation of cell culture media containing L-Glutamine that may be stored at 4°C for extended periods. The dipeptide is metabolized within the cells to yield L-Glutamine plus a second amino acid.

Uses:

Supplements, such as antibiotics, should be added as sterile supplements to the medium.

Storage conditions and shelf-life of supplemented product will be affected by the nature of the supplements.