

FETAL BOVINE SERUM ULTRA-LOW ENDOTOXIN

We are the preferred choice for all research, academic and industrial applications, which requires lowest endotoxin levels.

- | Lowest endotoxin level : < 0,1 EU/mL
- | Triple 0,1 um Filtered
- | Virus & Mycoplasma Tested
- | Ideal for sensitive Cells & Applications

All sera are to tested to determine and quantify endotoxin levels. Serox conducts a chromogenic quantitative test, method D of the European Pharmacopeia.

